

## GUYON KIDNEY PEDICLE CLAMP

curved; serrated

24.0cm (9.6")



**SKU:** 040-509-24

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Kidney Pedicle Clamp Forceps</a>
<b>Aesculap</b>	BH839
<b>KLS Martin</b>	13-531-24
<b>AUTHORS</b>	<a href="#">GUYON</a>
<b>Medicon</b>	51.56.24
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Kidney Pedicle Clamp Forceps</a>
<b>Aesculap</b>	BH839
<b>KLS Martin</b>	13-531-24

## AUTHORS

[GUYON](#)

## Medicon

51.56.24

## SPECIFICATIONS

[curved](#), [serrated](#)