

GUYON KIDNEY PEDICLE CLAMP

curved; serrated

23.0cm (9.2")




SKU: 040-513-23

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Kidney Pedicle Clamp Forceps
AUTHORS	GUYON
Medicon	51.56.15
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Kidney Pedicle Clamp Forceps
AUTHORS	GUYON
Medicon	51.56.15
SPECIFICATIONS	curved , serrated