

## GUYON KIDNEY PEDICLE CLAMP

curved; serrated

20.0cm (8")



**SKU:** 040-509-20

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Kidney Pedicle Clamp Forceps</a>
<b>Aesculap</b>	BH835
<b>KLS Martin</b>	13-531-20
<b>AUTHORS</b>	<a href="#">GUYON</a>
<b>Medicon</b>	51.56.20
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Kidney Pedicle Clamp Forceps</a>
<b>Aesculap</b>	BH835
<b>KLS Martin</b>	13-531-20

**AUTHORS**

[GUYON](#)

**Medicon**

51.56.20

**SPECIFICATIONS**

[curved](#), [serrated](#)