

## GUYON PEDICLE CLAMP

atraumatic; curved

24.0cm (9.6")




**SKU:** 040-519-24

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>                         |
| <b>TYPE</b>           | <a href="#">Kidney Pedicle Clamp Forceps</a> |
| <b>Aesculap</b>       | EF001  |
| <b>KLS Martin</b>     | 13-529-23                                    |
| <b>AUTHORS</b>        | <a href="#">GUYON</a>                        |
| <b>Medicon</b>        | 51.56.25                                     |
| <b>SPECIFICATIONS</b> | <a href="#">curved</a>                       |



### ZUSÄTZLICHE INFORMATIONEN

|                   |  |
|-------------------|--|
| <b>UNIT</b>       | <a href="#">pcs.</a>                         |
| <b>TYPE</b>       | <a href="#">Kidney Pedicle Clamp Forceps</a> |
| <b>Aesculap</b>   | EF001  |
| <b>KLS Martin</b> | 13-529-23                                    |

**AUTHORS**

[GUYON](#)

**Medicon**

51.56.25

**SPECIFICATIONS**

[curved](#)