

GYNECOLOGY/UROLOGY TELESCOPE HD

autoclavable

Ø 2.0mm; 30°; 300 mm

SKU: 601-132-30

Categories: [Endoscopy](#), [Equipment](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Gynecology](#)