

## GYNECOLOGY/UROLOGY TELESCOPE HD

autoclavable

Ø 2.0mm; 30°; 300 mm

**SKU:** 601-132-30

**Categories:** [Endoscopy](#), [Equipment](#), [Gynecology](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Gynecology](#)



30°

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Telescopes](#)

**SCOPE** [Gynecology](#)