

## HALSEY NEEDLE HOLDER

straight; serrated

13.0cm (5.2")

**SKU:** 050-040-13

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM205
<b>AUTHORS</b>	<a href="#">HALSEY</a>
<b>Medicon</b>	10.12.13
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM205
<b>AUTHORS</b>	<a href="#">HALSEY</a>
<b>Medicon</b>	10.12.13

## SPECIFICATIONS

[serrated](#), [straight](#)