

## HALSTED HAEMOSTATIC FORCEPS

straight; serrated

21.0cm (8.4")



**SKU:** 040-050-21

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                |
| <b>TYPE</b>           | <a href="#">Artery Forceps</a>                      |
| <b>Aesculap</b>       | BH210   |
| <b>KLS Martin</b>     | 13-316-21   |
| <b>AUTHORS</b>        | <a href="#">HALSTED</a>                             |
| <b>Medicon</b>        | 16.40.21  |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                   |                                |
|-------------------|--------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>           |
| <b>TYPE</b>       | <a href="#">Artery Forceps</a> |
| <b>Aesculap</b>   | BH210                          |
| <b>KLS Martin</b> | 13-316-21                      |

## AUTHORS

[HALSTED](#)

## Medicon

16.40.21

## SPECIFICATIONS

[serrated](#), [straight](#)