

HALSTED HAEMOSTATIC FORCEPS

straight; serrated

21.0cm (8.4")

SKU: 040-050-21

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH210
KLS Martin	13-316-21
AUTHORS	HALSTED
Medicon	16.40.21
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH210
KLS Martin	13-316-21

AUTHORS

[HALSTED](#)

Medicon

16.40.21

SPECIFICATIONS

[serrated](#), [straight](#)