

HALSTED HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

18.0cm (7.2")



SKU: 040-070-18

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH212
KLS Martin	13-318-18
AUTHORS	HALSTED
Medicon	16.42.18
SPECIFICATIONS	serrated , straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH212
KLS Martin	13-318-18

AUTHORS

[HALSTED](#)

Medicon

16.42.18

SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)