

## HALSTED HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

21.0cm (8.4")

**SKU:** 040-070-21

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH214
<b>KLS Martin</b>	13-318-21
<b>AUTHORS</b>	<a href="#">HALSTED</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH214
<b>KLS Martin</b>	13-318-21
<b>AUTHORS</b>	<a href="#">HALSTED</a>

## SPECIFICATIONS

[serrated](#), [straight](#), [toothed](#)