

HALSTED HAEMOSTATIC FORCEPS

straight; serrated

18.0cm (7.2")



SKU: 040-050-18

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH202
KLS Martin	13-316-18
AUTHORS	HALSTED
Medicon	16.40.18
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH202
KLS Martin	13-316-18

AUTHORS

[HALSTED](#)

Medicon

16.40.18

SPECIFICATIONS

[serrated](#), [straight](#)