

## HALSTED HAEMOSTATIC FORCEPS

straight; serrated

18.0cm (7.2")



**SKU:** 040-050-18

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH202
<b>KLS Martin</b>	13-316-18
<b>AUTHORS</b>	<a href="#">HALSTED</a>
<b>Medicon</b>	16.40.18
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH202
<b>KLS Martin</b>	13-316-18

## AUTHORS

[HALSTED](#)

## Medicon

16.40.18

## SPECIFICATIONS

[serrated](#), [straight](#)