


## HANDLE FOR POLYPECTOMY LOOPS

solid,stroke 200mm, rotating thumb ring

HF-connector 3mm



**SKU:** PK-4PWPk3000200

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Handles, Accessories](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Handles, Accessories](#)