

## HANK UTERINE DILATOR

double ended; fig. 19/20 (Ø 9.5 / 10.0mm)

**SKU:** 180-360-10

**Categories:** [Gynecology](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Uterine Dilators</a>
<b>KLS Martin</b>	32-460-20
<b>AUTHORS</b>	<a href="#">HANK</a>
<b>Medicon</b>	52.37.19
<b>SPECIFICATIONS</b>	<a href="#">double ended</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Uterine Dilators</a>
<b>KLS Martin</b>	32-460-20
<b>AUTHORS</b>	<a href="#">HANK</a>
<b>Medicon</b>	52.37.19
<b>SPECIFICATIONS</b>	<a href="#">double ended</a>