

## HARTMANN-HOFFMANN EAR FORCEPS

cup jaw Ø 2mm; trough-cutting

shaft 8.5cm (3.4")



**SKU:** 111-298-20

**Categories:** [Open Surgery](#), [Otology](#)

### ZUSÄTZLICHE INFORMATIONEN

|                   |   |
|-------------------|---|
| <b>UNIT</b>       | <a href="#">pcs.</a>                                |
| <b>TYPE</b>       | <a href="#">Ear Forceps</a>                         |
| <b>KLS Martin</b> | 36-456-08   |
| <b>AUTHORS</b>    | <a href="#">HARTMANN</a> , <a href="#">HOFFMANN</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                   |   |
|-------------------|---|
| <b>UNIT</b>       | <a href="#">pcs.</a>                                |
| <b>TYPE</b>       | <a href="#">Ear Forceps</a>                         |
| <b>KLS Martin</b> | 36-456-08   |
| <b>AUTHORS</b>    | <a href="#">HARTMANN</a> , <a href="#">HOFFMANN</a> |