

## HEANEY NEEDLE HOLDER

curved; serrated

21.0cm (8.4")




**SKU:** 050-141-21

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM259
<b>KLS Martin</b>	20-403-21
<b>AUTHORS</b>	<a href="#">HEANEY</a>
<b>Medicon</b>	10.31.21
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM259
<b>KLS Martin</b>	20-403-21

**AUTHORS**

[HEANEY](#)

**Medicon**

10.31.21

**SPECIFICATIONS**

[curved](#), [serrated](#)