

HEGAR-BAUMGARTNER NEEDLE HOLDER

straight; serrated

14.0cm (5.6")



SKU: 050-048-14

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM215
KLS Martin	20-340-14
AUTHORS	BAUMGARTNER , HEGAR
Medicon	10.14.14
SPECIFICATIONS	serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM215
KLS Martin	20-340-14

AUTHORS

[BAUMGARTNER](#), [HEGAR](#)

Medicon

10.14.14

SPECIFICATIONS

[serrated](#), [straight](#)