

## HEGAR UTERINE DILATOR

with sloped handle; Ø 13.0mm

18.5cm (7.25")

**SKU:** 180-350-13

**Categories:** [Gynecology](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Uterine Dilators</a>
<b>Aesculap</b>	EM113
<b>KLS Martin</b>	32-400-13
<b>AUTHORS</b>	<a href="#">HEGAR</a>
<b>Medicon</b>	52.30.13
<b>SPECIFICATIONS</b>	<a href="#">sloped handle</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Uterine Dilators</a>
<b>Aesculap</b>	EM113
<b>KLS Martin</b>	32-400-13

## AUTHORS

[HEGAR](#)

## Medicon

52.30.13

## SPECIFICATIONS

[sloped handle](#)