

HEGAR UTERINE DILATOR

with sloped handle; Ø 13.0mm

18.5cm (7.25")

SKU: 180-350-13

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|----------------------------------|
| UNIT | pcs. |
| TYPE | Uterine Dilators |
| Aesculap | EM113 |
| KLS Martin | 32-400-13 |
| AUTHORS | HEGAR |
| Medicon | 52.30.13 |
| SPECIFICATIONS | sloped handle |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|----------------------------------|
| UNIT | pcs. |
| TYPE | Uterine Dilators |
| Aesculap | EM113 |
| KLS Martin | 32-400-13 |

AUTHORS

[HEGAR](#)

Medicon

52.30.13

SPECIFICATIONS

[sloped handle](#)