

HEGAR UTERINE DILATOR

with sloped handle; Ø 14.0mm

18.5cm (7.25")

SKU: 180-350-14

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM114
KLS Martin	32-400-14
AUTHORS	HEGAR
Medicon	52.30.14
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM114
KLS Martin	32-400-14

AUTHORS

[HEGAR](#)

Medicon

52.30.14

SPECIFICATIONS

[sloped handle](#)