

HEGAR UTERINE DILATOR

with sloped handle; Ø 15.0mm

18.5cm (7.25")

SKU: 180-350-15

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM115
KLS Martin	32-400-15
AUTHORS	HEGAR
Medicon	52.30.15
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM115
KLS Martin	32-400-15

AUTHORS

[HEGAR](#)

Medicon

52.30.15

SPECIFICATIONS

[sloped handle](#)