

HEGAR UTERINE DILATOR

with sloped handle; Ø 17.0mm

18.5cm (7.25")

SKU: 180-350-17

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM117
KLS Martin	32-400-17
AUTHORS	HEGAR
Medicon	52.30.17
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM117
KLS Martin	32-400-17

AUTHORS

[HEGAR](#)

Medicon

52.30.17

SPECIFICATIONS

[sloped handle](#)