

HEGAR UTERINE DILATOR

with sloped handle; Ø 24.0mm

18.5cm (7.25")

SKU: 180-350-24

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM124
KLS Martin	32-400-24
AUTHORS	HEGAR
Medicon	52.30.24
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM124
KLS Martin	32-400-24

AUTHORS

[HEGAR](#)

Medicon

52.30.24

SPECIFICATIONS

[sloped handle](#)