

HEGAR UTERINE DILATOR

with sloped handle; Ø 27.0mm

18.5cm (7.25")

SKU: 180-350-27

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-400-27
AUTHORS	HEGAR
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-400-27
AUTHORS	HEGAR
SPECIFICATIONS	sloped handle