

HEGAR UTERINE DILATOR

with sloped handle; Ø 28.0mm

18.5cm (7.25")

SKU: 180-350-28

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-400-28
AUTHORS	HEGAR
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-400-28
AUTHORS	HEGAR
SPECIFICATIONS	sloped handle