

## HEGAR UTERINE DILATOR

with sloped handle; Ø 3.5mm

18.5cm (7.25")

**SKU:** 180-351-03

**Categories:** [Gynecology](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Uterine Dilators</a>
<b>Aesculap</b>	EM133
<b>KLS Martin</b>	32-401-03
<b>AUTHORS</b>	<a href="#">HEGAR</a>
<b>Medicon</b>	52.31.03
<b>SPECIFICATIONS</b>	<a href="#">sloped handle</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Uterine Dilators</a>
<b>Aesculap</b>	EM133
<b>KLS Martin</b>	32-401-03

## AUTHORS

[HEGAR](#)

## Medicon

52.31.03

## SPECIFICATIONS

[sloped handle](#)