

HEGAR UTERINE DILATOR

with sloped handle; Ø 9.5mm

18.5cm (7.25")

SKU: 180-351-09

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM139
KLS Martin	32-401-09
AUTHORS	HEGAR
Medicon	52.31.09
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM139
KLS Martin	32-401-09

AUTHORS

[HEGAR](#)

Medicon

52.31.09

SPECIFICATIONS

[sloped handle](#)