

HEGAR UTERINE DILATOR

with sloped handle; Ø 15.5mm

18.5cm (7.25")

SKU: 180-351-15

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|----------------------------------|
| UNIT | pcs. |
| TYPE | Uterine Dilators |
| Aesculap | EM145 |
| KLS Martin | 32-401-15 |
| AUTHORS | HEGAR |
| Medicon | 52.31.15 |
| SPECIFICATIONS | sloped handle |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|----------------------------------|
| UNIT | pcs. |
| TYPE | Uterine Dilators |
| Aesculap | EM145 |
| KLS Martin | 32-401-15 |

AUTHORS

[HEGAR](#)

Medicon

52.31.15

SPECIFICATIONS

[sloped handle](#)