

HEGAR UTERINE DILATOR

with sloped handle; Ø 18.5mm

18.5cm (7.25")

SKU: 180-351-18

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-401-18
AUTHORS	HEGAR
Medicon	52.31.18
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-401-18
AUTHORS	HEGAR
Medicon	52.31.18

SPECIFICATIONS

[sloped handle](#)