

HEGAR UTERINE DILATOR

with sloped handle; Ø 19.5mm

18.5cm (7.25")

SKU: 180-351-19

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-401-19
AUTHORS	HEGAR
Medicon	52.31.19
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-401-19
AUTHORS	HEGAR
Medicon	52.31.19

SPECIFICATIONS

[sloped handle](#)