

HEGAR UTERINE DILATOR

with sloped handle; Ø 5.0mm

18.5cm (7.25")

SKU: 180-350-05

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM105
KLS Martin	32-400-05
AUTHORS	HEGAR
Medicon	52.30.05
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM105
KLS Martin	32-400-05

AUTHORS

[HEGAR](#)

Medicon

52.30.05

SPECIFICATIONS

[sloped handle](#)