

HEGAR UTERINE DILATOR

double ended; Ø 5.0 / 6.0mm

20.0cm (8")

SKU: 180-354-05

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM233
KLS Martin	32-430-06
AUTHORS	HEGAR
Medicon	52.35.05
SPECIFICATIONS	double ended

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM233
KLS Martin	32-430-06

AUTHORS

[HEGAR](#)

Medicon

52.35.05

SPECIFICATIONS

[double ended](#)