

HEGAR UTERINE DILATOR

double ended; Ø 11.0 / 12.0mm

20.0cm (8")

SKU: 180-354-11

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM236
KLS Martin	32-430-12
AUTHORS	HEGAR
Medicon	52.35.11
SPECIFICATIONS	double ended

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM236
KLS Martin	32-430-12

AUTHORS

[HEGAR](#)

Medicon

52.35.11

SPECIFICATIONS

[double ended](#)