

## HEGAR UTERINE DILATOR

double ended; Ø 11.0 / 12.0mm

20.0cm (8")

**SKU:** 180-354-11

**Categories:** [Gynecology](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Uterine Dilators</a>
<b>Aesculap</b>	EM236
<b>KLS Martin</b>	32-430-12
<b>AUTHORS</b>	<a href="#">HEGAR</a>
<b>Medicon</b>	52.35.11
<b>SPECIFICATIONS</b>	<a href="#">double ended</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Uterine Dilators</a>
<b>Aesculap</b>	EM236
<b>KLS Martin</b>	32-430-12

**AUTHORS**

[HEGAR](#)

**Medicon**

52.35.11

**SPECIFICATIONS**

[double ended](#)