

## HEGAR UTERINE DILATOR

double ended; Ø 13.0 / 14.0mm

20.0cm (8")

**SKU:** 180-354-13

**Categories:** [Gynecology](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

|                       |                                  |
|-----------------------|----------------------------------|
| <b>UNIT</b>           | <a href="#">pcs.</a>             |
| <b>TYPE</b>           | <a href="#">Uterine Dilators</a> |
| <b>Aesculap</b>       | EM237                            |
| <b>KLS Martin</b>     | 32-430-14                        |
| <b>AUTHORS</b>        | <a href="#">HEGAR</a>            |
| <b>Medicon</b>        | 52.35.13                         |
| <b>SPECIFICATIONS</b> | <a href="#">double ended</a>     |

### ZUSÄTZLICHE INFORMATIONEN

|                   |                                  |
|-------------------|----------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>             |
| <b>TYPE</b>       | <a href="#">Uterine Dilators</a> |
| <b>Aesculap</b>   | EM237                            |
| <b>KLS Martin</b> | 32-430-14                        |

**AUTHORS**

[HEGAR](#)

**Medicon**

52.35.13

**SPECIFICATIONS**

[double ended](#)