

HEGAR UTERINE DILATOR

double ended; Ø 15.0 / 16.0mm

20.0cm (8")

SKU: 180-354-15

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM238
KLS Martin	32-430-16
AUTHORS	HEGAR
Medicon	52.35.15
SPECIFICATIONS	double ended

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM238
KLS Martin	32-430-16

AUTHORS

[HEGAR](#)

Medicon

52.35.15

SPECIFICATIONS

[double ended](#)