

HEGAR UTERINE DILATOR

double ended; Ø 17.0 / 18.0mm

20.0cm (8")

SKU: 180-354-17

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM239
KLS Martin	32-430-18
AUTHORS	HEGAR
Medicon	52.35.17
SPECIFICATIONS	double ended

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM239
KLS Martin	32-430-18

AUTHORS

[HEGAR](#)

Medicon

52.35.17

SPECIFICATIONS

[double ended](#)