

HEGAR UTERINE DILATOR

double ended; Ø 23.0 / 24.0mm

20.0cm (8")

SKU: 180-354-23

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-430-24
AUTHORS	HEGAR
Medicon	52.35.23
SPECIFICATIONS	double ended

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-430-24
AUTHORS	HEGAR
Medicon	52.35.23

SPECIFICATIONS

[double ended](#)