

HEGAR UTERINE DILATOR

double ended; Ø 25.0 / 26.0mm

20.0cm (8")

SKU: 180-354-25

Categories: [Gynecology](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-430-26
AUTHORS	HEGAR
Medicon	52.35.25
SPECIFICATIONS	double ended

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-430-26
AUTHORS	HEGAR
Medicon	52.35.25

SPECIFICATIONS

[double ended](#)