

**HEGAR UTERINE DILATOR**

with sloped handle; Ø 1.0mm  
 18.5cm (7.25")

**SKU:** 180-350-01  
**Categories:** [Gynecology](#), [Open Surgery](#)



**ZUSÄTZLICHE INFORMATIONEN**

|                       |                                  |
|-----------------------|----------------------------------|
| <b>UNIT</b>           | <a href="#">pcs.</a>             |
| <b>TYPE</b>           | <a href="#">Uterine Dilators</a> |
| <b>Aesculap</b>       | EM1001                           |
| <b>KLS Martin</b>     | 32-400-01                        |
| <b>AUTHORS</b>        | <a href="#">HEGAR</a>            |
| <b>Medicon</b>        | 52.30.01                         |
| <b>SPECIFICATIONS</b> | <a href="#">sloped handle</a>    |

**ZUSÄTZLICHE INFORMATIONEN**

|                   |                                  |
|-------------------|----------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>             |
| <b>TYPE</b>       | <a href="#">Uterine Dilators</a> |
| <b>Aesculap</b>   | EM1001                           |
| <b>KLS Martin</b> | 32-400-01                        |

## AUTHORS

[HEGAR](#)

## Medicon

52.30.01

## SPECIFICATIONS

[sloped handle](#)