

HEGAR UTERINE DILATORS SET

with 26 pieces; Ø 1-26mm; in metal case



SKU: 180-350-99

Categories: [Gynecology](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Uterine Dilators](#)

Aesculap EM192

KLS Martin 32-400-90

AUTHORS [HEGAR](#)

Medicon 52.30.45

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Uterine Dilators](#)

Aesculap EM192

KLS Martin 32-400-90

AUTHORS [HEGAR](#)

Medicon 52.30.45