

HEISS HAEMOSTATIC FORCEPS

slightly curved; serrated

20.0cm (8")




SKU: 040-391-20

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Artery Forceps , Bronchus Forceps |
| Aesculap | BH207 |
| KLS Martin | 13-423-20 |
| AUTHORS | HEISS |
| Medicon | 16.51.20 |
| SPECIFICATIONS | curved , serrated |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---|
| UNIT | pcs. |
| TYPE | Artery Forceps , Bronchus Forceps |
| Aesculap | BH207 |
| KLS Martin | 13-423-20 |

AUTHORS

[HEISS](#)

Medicon

16.51.20

SPECIFICATIONS

[curved](#), [serrated](#)