

## HEISS HAEMOSTATIC FORCEPS

straight; serrated

20.0cm (8")



**SKU:** 040-390-20

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH206
<b>KLS Martin</b>	13-422-20
<b>AUTHORS</b>	<a href="#">HEISS</a>
<b>Medicon</b>	16.50.20
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH206
<b>KLS Martin</b>	13-422-20

## AUTHORS

[HEISS](#)

## Medicon

16.50.20

## SPECIFICATIONS

[serrated](#), [straight](#)