

HENKE TONSIL DISSECTOR

double ended; 16.5mm / 5.0mm

24.5cm (9.8")

SKU: 132-456-23

Categories: [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tonsil Dissectors
Aesculap	OM669
KLS Martin	39-291-23
AUTHORS	HENKE



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tonsil Dissectors
Aesculap	OM669
KLS Martin	39-291-23
AUTHORS	HENKE