

HENKE TONSIL DISSECTOR

double ended; 11.5mm / 5.0mm

24.5cm (9.8")

SKU: 132-455-23

Categories: [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tonsil Dissectors
Aesculap	OM668
KLS Martin	39-289-23
AUTHORS	HENKE



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tonsil Dissectors
Aesculap	OM668
KLS Martin	39-289-23
AUTHORS	HENKE