

## HENKE TONSIL DISSECTOR

double ended; 11.5mm / 5.0mm

24.5cm (9.8")

**SKU:** 132-455-23

**Categories:** [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tonsil Dissectors</a>
<b>Aesculap</b>	OM668
<b>KLS Martin</b>	39-289-23
<b>AUTHORS</b>	<a href="#">HENKE</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tonsil Dissectors</a>
<b>Aesculap</b>	OM668
<b>KLS Martin</b>	39-289-23
<b>AUTHORS</b>	<a href="#">HENKE</a>