

## HENLEY SELFRETAINING RETRACTOR

blunt; 3X4 teeth; complete with blades

16.5cm (6.6")




**SKU:** 070-301-00

**Categories:** [Open Surgery](#), [Retractors](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                            |
| <b>TYPE</b>           | <a href="#">Selfretaining Retractors</a>        |
| <b>KLS Martin</b>     | 15-768-00                                       |
| <b>AUTHORS</b>        | <a href="#">HENLEY</a>                          |
| <b>SPECIFICATIONS</b> | <a href="#">blunt</a> , <a href="#">toothed</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                            |
| <b>TYPE</b>           | <a href="#">Selfretaining Retractors</a>        |
| <b>KLS Martin</b>     | 15-768-00                                       |
| <b>AUTHORS</b>        | <a href="#">HENLEY</a>                          |
| <b>SPECIFICATIONS</b> | <a href="#">blunt</a> , <a href="#">toothed</a> |