

HENLEY SELFRETAINING RETRACTOR

blunt; 3X4 teeth; complete with blades

16.5cm (6.6")




SKU: 070-301-00

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-768-00
AUTHORS	HENLEY
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-768-00
AUTHORS	HENLEY
SPECIFICATIONS	blunt , toothed