

## HENLEY SELFRETAINING RETRACTOR ONLY

blunt; 3X4 teeth

16.0cm (6.4")




**SKU:** 070-301-16

**Categories:** [Open Surgery](#), [Retractors](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>AUTHORS</b>	<a href="#">HENLEY</a>
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>AUTHORS</b>	<a href="#">HENLEY</a>
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">toothed</a>